



TPR100 2X

PRODUCT CODE:

D173, D179,
BE175

DESCRIPTION:

Specially formulated hard water stable presoak. Hard water stable up to 10 grains hardness at 200:1 dilution.

PRODUCT APPLICATION:

Presoak

DILUTION RECOMMENDATIONS:

Presoak

2x
120-240:1

2x

Step 1: Fill graduated cylinder to **10 ml.** mark with solution.

Step 2: Add **4** drops Indicator **P.** Solution will turn pink.

Step 3: Add Titrate **#1** dropwise. Invert tube to mix after each drop.

Step 4: Continue adding drops until solution turns colorless.

Dilution Range: 3001 – 50:1

Indicator: P

Titrate #: 1



| <u># DROPS</u> | <u>OZ. PER GAL.</u> | <u>RATIO</u> |
|----------------|---------------------|--------------|
| 3 | 0.43 | 300 |
| 4 | 0.59 | 217 |
| 5 | 0.76 | 170 |
| 6 | 0.92 | 139 |
| 7 | 1.08 | 118 |
| 8 | 1.25 | 103 |
| 9 | 1.41 | 91 |
| 10 | 1.58 | 81 |
| 11 | 1.74 | 74 |
| 12 | 1.90 | 67 |
| 13 | 2.07 | 62 |
| 14 | 2.23 | 57 |
| 15 | 2.40 | 53 |
| 16 | 2.56 | 50 |